



Sieve Analysis

Different sieving methods for a variety of applications

AUTHOR

Mike Lucka



Product Manager

Retsch[®]
Solutions in Milling & Sieving

Retsch GmbH
Retsch-Allee 1-5
42781 Haan, Germany

Phone: +49 (0)2104/2333-100
E-Mail: m.lucka@retsch.com

www.retsch.com

The solubility of drugs, the taste of chocolate or the solidity of concrete – these are only a few examples of physical or chemical material properties which are influenced by particle size. Hence, the determination and knowledge of the particle size distribution is an essential part of the quality control process for industrial products. From incoming and production control to research and development sieve analyses are used to determine a number of parameters or simply the particle size. Easy handling, low investment cost and high accuracy make sieve analysis one of the most frequently used procedures for measuring the particle size. This white paper gives an overview of the different sieving techniques and describes the necessary steps to ensure reliable results.


Sieving is carried out to separate a sample according to its particle sizes by submitting it to mechanical force. The direction, intensity and type of force depends on the chosen sieving method. The sample is either moved in horizontal or vertical direction. In tap sieve shakers both movements are superimposed. Air jet sieving is a special case, here the sample is dispersed by an air jet blown out of a rotating nozzle.

Which factors need to be established to select a suitable sieving method?

- Particle size

If the required measuring range lies between approximately 40 microns and 125 mm, classical dry sieving is the method of choice. The range may be extended to 20 microns by wet sieving and to 10 microns by air jet sieving (figure 1).

Fig. 1:
Measuring range of air jet, wet and dry sieving



● **Sample properties**

It should be taken into account whether the particles form agglomerates, which density the material has or if it tends to be electrostatically charged.

● **Standards**

Sieving methods are described in the DIN standard 66165. If industry-specific test procedures or standards exist these also determine the choice of method.

● **Number of fractions**

Are several fractions required? Or is it sufficient to know which percentage of the sample is smaller or bigger than a defined particle size? The method for obtaining the latter information is called a **sieve cut** because the sample is simply separated into two fractions.

Different Sieving Methods

Once these questions are satisfactorily answered, a suitable sieving method can be selected. This white paper describes the different sieving techniques.

Sieving with controlled acceleration

The Sieve Shakers AS 200 control, AS 300 control and AS 450 control are activated in their natural frequency. This means, the sieving frequency changes with the load of the instrument. It depends on the weight of the sieve stack and the sample quantity. In order to ensure the reproducibility of the results even in short-time sieving procedures, the default setting of the vibration height can be switched to sieve acceleration (sieving with equal acceleration).

The RETSCH Sieve Shakers AS 200 control, AS 300 control and AS 450 control are the only sieve shakers to feature the possibility of eliminating influences of error by differing sieving frequencies by automatic amplitude adjustment.

Vibrational sieving (figure 2): Vibrational sieving submits the sample to **three dimensional movements**. A vertical throwing motion is superimposed by a circular movement. This mechanism causes the particles to be evenly distributed over the entire sieving surface and to be thrown into the air where they ideally change their orientation in a way that allows them to be compared to the sieve apertures in all possible dimensions. Vibratory sieve shakers like the AS 200 series (figure 3), AS 300 control and AS 450 basic and control operate according to this principle. All RETSCH vibratory sieve shakers are suitable for **dry as well as wet sieving**. The "control" models can be calibrated and permit software-based evaluation of the sieving process. Moreover, they provide reproducible, globally comparable results thanks to **sieving with controlled acceleration** (see box on the left).



Fig. 2 (left):
Principle of vibratory sieving

Fig. 3 (right):
Vibratory sieve shakers of the AS 200 series

	AS 200 basic	AS 200 digit cA	AS 200 control	AS 300 control	AS 450 basic	AS 450 control
Measuring range	20 µm - 25 mm			20 µm - 40 mm	25 µm - 125 mm	25 µm - 125 mm
Sieve diameter	100 / 150 / 200 / 203 (8") mm			100 / 150 / 200 / 203 (8") / 305 (12") / 315 mm	400 / 450 mm	400 / 450 mm
Amplitude in mm	1 - 100%	0.2 - 3.0	0.2 - 3.0	0.2 - 2.2	0.2 - 2	0 - 2.2
Maximum load	3 kg			6 kg	15 kg	25 kg
Recalibration and sieve acceleration	-	-	yes	Yes	-	yes

Table 1:
A suitable vibratory sieve shaker for every requirement

Wet sieving: A special case of vibratory sieving is wet sieving. Agglomerates, electrostatic charging or a high degree of fineness can all make the sieving process difficult and in such cases wet sieving may be called for. This involves washing the sample through the sieve stack in a suitable medium which is usually water.

Horizontal sieving (figure 4): Here the sample is subjected to a **circular horizontal movement**. This two-dimensional motion does not cause the particles to change their original orientation. This method is particularly suitable for **longish, disk-shaped or fibrous samples** (sieving with circular motion according to DIN 53 477). For this application RETSCH offers the AS 400 control with an electronically controlled speed of 50 to 300 rpm. Sieve diameters of 100 mm / 150 mm / 200 mm / 203 mm (8") / 305 (12") mm / 315 mm / 400 mm may be used to separate bulk samples of up to 15 kg in one run. The measuring range lies between 45 µm and 63 mm. This sieve shaker can be recalibrated and controlled with RETSCH's evaluation software EasySieve[®].



Fig. 4 (left):
Principle of horizontal sieving



Fig. 5 (right):
Horizontal sieve shaker
AS 400 control

Tap sieving (figure 6): The use of tap sieve shakers is stipulated in a number of standards. In the tap sieving process, a circular horizontal movement is superimposed by a vertical tapping motion as, for example, in the AS 200 tap (speed 280 rpm, taps 150 per min). The AS 200 tap (fig. 7) accepts up to 7 sieves with diameters of 200 mm or 203 mm (8") and a maximum sample load of 3 kg. This sieve shaker can also be controlled with the evaluation software EasySieve[®].



Fig. 6 (left):
Principle of tap sieving

Fig. 7 (right):
Tap sieve shaker AS 200 tap



Air jet sieving (figure 8): To obtain a **sieve cut** by air jet sieving **only one sieve** is used instead of a stack of sieves, and the sieve itself is not put into motion. An industrial vacuum cleaner generates low pressure inside the sieve chamber. The sucked-in air escapes with high speed from the rotating slit nozzle below the sieve and disperses the particles which can then be compared to the sieve apertures. When the particles hit the sieve lid they are not only redirected but also **deagglomerated**. The particles which are small enough are then transported through the sieve mesh and sucked in. It is also possible to collect them in a **cyclone**. RETSCH offers the AS 200 jet for air jet sieving (figure 9). Features include the innovative **Open Mesh Function** to reduce the number of near-mesh particles and automatic vacuum regulation (accessory). Sieve diameters of 200 mm and 203 mm (8") can be used with the AS 200 jet which covers a measuring range from 10 µm to 4 mm. The machine may be controlled with EasySieve[®] evaluation software.




Fig. 8 (left):
Principle of air jet sieving

Fig. 9 (right):
Air jet sieving machine
AS 200 jet



Sieve Analysis Step by Step

A complete sieving process consists of the following steps which should be carried out precisely and carefully.

- a) Sampling
- b) Sample division (if required)
- c) Selection of suitable test sieves
- d) Selection of sieving parameters
- e) Actual sieve analysis
- f) Recovery of sample material
- g) Data evaluation
- h) Cleaning and drying the test sieves

a) Sampling

Correct sampling is essential to obtain a representative laboratory sample, especially if the material is heterogeneous. Some important questions should be clarified in advance:

- Which quantity is required to ensure representativeness of the original sample
- From which part of the original material should the sample be taken

Some industry-specific standards provide guidelines and directions on the correct sampling process, for example DIN 51701 in the coal industry. It contains, for example, the formula

$$G \text{ [kg]} = 0.07 \text{ [kg/mm]} \times z \text{ [mm]}$$

which indicates how much sample "G" must be extracted from a bulk sample with maximum particle size "z" to obtain a representative quantity. Taking a coal sample with a maximum particle size of 5 cm as an example, the following calculation applies:

$$G \text{ [kg]} = 0.07 \text{ kg/mm} \times 50 \text{ mm}$$

$$G = 3.5 \text{ kg}$$

Consequently, the extracted sample amount should be at least 3.5 kg to make sure it is representative. For other materials than coal a different density should be taken into account.

b) Sample division

If a sample quantity is too large for analysis, it must be divided in a way that ensures the part sample represents the original material. Sampling may be carried out randomly or **manually** by coning and quartering. These methods, however, are **error-prone** and may cause standard deviations amounting to 10% and more just by imprecise division (figure 10, C + D). **Automated sample dividers** like the PT 100 (figure 11) or PT 200 from RETSCH offer a much more exact and convenient method. In the PT 100 the sample flows through a feed hopper directly into the openings of a dividing head. Even if the sample is coarse-grained, the deviations among the quantities of each sample vessel are negligible. The division process runs fully automatically. The dividing head – which can have 6, 8 or 10 outlets – executes a controlled constant 11 rotations per minute, independent of load and power frequency. This means practically that the sample flow is divided 1100 times per minute when a dividing head with 10 outlets is used. Thus the **highest possible degree of division accuracy** is statistically guaranteed.



Fig. 10:
Rotating sample divider PT 100




Fig. 11: Standard deviations of analysis results resulting from different sampling methods

c) Selection of suitable test sieves

The selection of the sieves depends on the **sample quantity** and on the approximate **particle size distribution**. To cover the complete particle size range of the sample the increments of the mesh apertures should follow a logarithmic series. This is frequently done with the help of the Renard series, for example R10/3. The standard DIN 66165 contains some guidelines on the selection of sieve diameters and mesh sizes. It also stipulates the maximum sample load permitted for different mesh sizes and also the maximum particle size (figure 12).

Correct loading of the sieves ensures reproducible results!

Calculation of sieve load

The oversize on a sieve with a mesh size of 1 mm, for example, should not be more than 20 cm³ per square decimeter. For a 200 mm sieve that equals 63 cm³ oversize, for a 400 mm sieve it is 252 cm³. The maximum batch should not exceed twice the amount of the oversize value, i.e. a 200 mm sieve with mesh size 1 mm should not be filled with more than 126 cm³ sample material. By multiplying these values with the bulk density, the corresponding masses can be obtained.

Mesh size	Max. batch	Max. permitted sieve oversize
25 µm	14 cm ³	7 cm ³
45 µm	20 cm ³	10 cm ³
63 µm	26 cm ³	13 cm ³
125 µm	38 cm ³	19 cm ³
250 µm	58 cm ³	29 cm ³
500 µm	88 cm ³	44 cm ³
1 mm	126 cm ³	63 cm ³
2 mm	220 cm ³	110 cm ³
4 mm	346 cm ³	173 cm ³
8 mm	566 cm ³	283 cm ³

Fig. 12: Examples for the maximum batch and permitted sieve oversize for 200 mm sieves (according to DIN 66165)

The maximum sieve oversize permitted on a sieve can be approximately calculated with the formula

$$R = 0.00178 \times D^2 \times w^{0.667} \times \rho$$

with R being the maximum permitted sieve oversize in [g]. D stands for the sieve diameter (200 mm / 300 mm etc.), w is the mesh size in [mm] and ρ describes the material density. In the standard 66165 and in the above example a density of 1 kg/dm³ was assumed. The standard also describes the maximum load per sieve which should not exceed the maximum sieve oversize D by more than twice the amount.

The maximum particle size permitted can be calculated by using this formula

$x_{max} = 10 \times w^{0.7}$ with x_{max} [mm] being the maximum particle size permitted, and w describing the nominal mesh size of each sieve (the value is used in the formula without unit). For a sample on a sieve with an aperture size of 250 microns the following calculation is used to determine the maximum particle size permitted:

$$x_{max} = 10 \times w^{0.7}$$

$$x_{max} = 10 \times 0.250^{0.7}$$

$$x_{max} = 3.79 \text{ (unit = mm)}$$

d) Selecting suitable sieving parameters

Usually, industry-specific standards or company-specific regulations contain information about the required sieving parameters. Should this not be the case, those parameters need to be ascertained by experiment. The sieving time and, when using vibratory sieve shakers also the amplitude, are the most important factors. Observing the maximum sieve load helps protecting the sieve from damages and ensures that every particle has the possibility to compare with the sieve mesh as frequently as possible and in every dimension.

Amplitude:

During a sieving process there is a continuous size comparison between sieve apertures and particles. If a particle is smaller than an aperture, it passes the mesh. If it does not pass the aperture it is thrown upwards with the next lifting of the sieve bottom, takes a different orientation – which is important with longish particles – and hits the mesh again for a new comparison. Each comparison is an opportunity for the particle to pass the sieve mesh. Hence, the objective of a sieving process is to **generate as many comparisons as possible in a given time**. This is ideally the case when within the period of a sieve bottom vibration exactly one comparison takes place and the particles which don't pass the sieve mesh are accelerated in a way that within the next period another comparison takes place. This state is called **statistical resonance** (blue curve in figure 13). If a particle is accelerated too strongly it has only very few occasions to compare with the sieve apertures (yellow curve) and the complete separation of the sample is delayed; if a particle is hardly accelerated or not at all, it cannot orientate freely for a new comparison (red curve). The acceleration rate of a particle depends on the set amplitude. **The optimum amplitude for a particular material and quantity needs to be ascertained empirically**. Amplitude here means the complete oscillation in the horizontal plane, below as well as above the idle position. An amplitude of 1.2 mm indicates a movement of + 0.6 mm and – 0.6 mm in relation to the idle position.




Fig. 13:
Movement of particles in a
vibratory sieve shaker in
relation to the sieve bottom

Blue curve: The particle is in statistical resonance with the sieve bottom; red curve: the particle is not sufficiently moved; yellow curve: the particle is thrown too high

Sieving time:

According to DIN 66165 the sieving process is considered as finished when, after one minute of sieving, less than 0.1% of the feed quantity passes the sieve. If the undersize is larger the sieving time needs to be extended. Experience has shown that vibratory sieving in 3 dimensions is particularly suitable for achieving a sharp separation of fractions in a very short time.

e) Carrying out the actual sieve analysis

After selection of the parameters the actual sieving process starts. The following steps have to be carried out in chronological order:

1. Put together a sieve stack with collecting pan
2. Select sieving aids, if required: for mesh sizes < 500 microns the use of cubes, chains or brushes is recommended to facilitate passage of the sample

3. Determine the empty weight of sieves and collecting pan: This can be done software-based or manually by using a balance. Suitable programs such as EasySieve[®] make weighing and evaluation much easier.
4. Place the sieve stack with increasing mesh size on the collecting pan (fig. 14)
5. Weigh the sample and put it on the top sieve (biggest mesh size); clamp sieve stack on the machine
6. Set the amplitude / speed and sieving time on the sieve shaker
7. Start the sieve shaker
8. When the sieving time has expired weigh each sieve and the collecting pan with the fraction on it
9. Determine the mass and percentage of each fraction
10. Evaluation



Fig. 14:
Stacking of sieves (example)

f) Recovery of sample

When the sieving process is finished the material is removed from the sieves. The recovery of the fractions is a decisive advantage of sieve analysis over most optical measurement methods. The fractions are not only analytical values but physically exist and can be used for further processes following sieve analysis.

g) Determination and evaluation of data

Data evaluation can be done manually or with the help of a software like EasySieve[®]. As shown in table 2 the percentages are calculated and displayed graphically (figure 15).

Sieve [µm]	Net weight [g]	Weight after sieving [g]	Difference [g]	Percentage p ₃ [%]	Cumulative distribution Q ₃ [%]
Pan	501	505.5	4.5	3	3
45	253	259	6	4	7
63	268	283	15	10	17
125	298	328	30	20	37
250	325	373	48	32	69
500	362	384.5	22.5	15	84
1,000	386	401	15	10	94
2,000	406	412	6	4	98
4,000	425	428	3	2	100
			= 150 g	= 100%	

Tabelle 2:
Example of a sieve analysis

The empty sieves are weighed before and, containing the respective fractions, after the sieving process. The difference corresponds to the weight of each fraction. If these are put into relation to the total weight of the sample the percentage of each fraction can be calculated. This method has the advantage that due to the absence of dimensions sieving is carried out independently of the density or mass of the sample material. The difference between weighed sample and sum of the single fractions is called sieve loss. If this is more than 1% of the feed quantity DIN 66165 stipulates that the sieve analysis needs to be repeated. The mass percentages of the fractions can be graphically displayed as histograms (figure 15). In the example the largest fractions can be found between 250 microns and 500 microns with 32%. By adding up the single fractions and through interpolation between the measuring points the cumulative distribution curve Q_3 is obtained (figure 16).

This curve helps to determine a variety of sample properties. Taking a look, for example, at the particle size 250 microns the corresponding value of 36% can be read off the y-axis. This means that 36% of the total sample are smaller than 250 microns (figure 17). To determine the median Q_3 (50) of this distribution, the corresponding particle size (330 microns) is read off the x-axis. This means that 50% of the sample mass are smaller than or equal 330 microns. In this way various $x(Q_3)$ or $Q_3(x)$ values can be determined.

A **software** like **EasySieve®** allows for **quick and error-free evaluation** of various parameters.

Fig. 15 (left):
Histogram with percentages of the fractions

Fig. 16 (right):
Histogram with cumulative distribution curve of the fractions






Fig. 17:
Cumulative distribution curve with exemplary percentages



h) Cleaning and drying of test sieves

Test sieves are measuring instruments which should be treated with care before, during and after sieving. By no means should the sample be forced through the sieve mesh during the sieving process. Even a light brushing of the material – particularly through very fine fabric – may lead to changes of the mesh and damage the sieve wire gauze. Near-mesh particles trapped in

the sieve mesh are easily removed by turning the sieve upside down and tapping it lightly on a table. For particles which cannot be removed this way a fine hair brush may be used to lightly brush the underside of the sieve mesh.

Coarser fabrics with mesh sizes larger than 500 microns can be effectively cleaned dry or wet with a hand brush with plastic bristles. Sieves with a mesh size below 500 microns should generally be cleaned in an **ultrasonic bath**. Drying cabinets of various sizes can be used for drying test sieves which should be placed vertically inside the cabinet. It is recommended not to exceed a temperature of 80 °C. With higher temperatures especially the fine metal wire mesh could become warped; as a result, the tension of the fabric inside the sieve frame is reduced which makes the sieve less efficient. **RETSCH's Fluid Bed Dryer TG 200** is particularly effective in drying test sieves with a diameter of 200/203 mm. Before cleaning or drying the sieves, the rubber or plastic seal rings have to be removed.

If deviations from the uniformity of the mesh are observed, the sieve is no longer suitable for quality control purposes (DIN ISO 9000 ff) and must be replaced.



Fig. 18:
Fluid Bed Dryer TG 200

Conclusion

Sieve analysis continues to be one of the best proven methods for particle size analysis of bulk materials and is mentioned in various standards. It is easy to perform and provides precise and reproducible results while pre-serving the physical size fractions. Moreover, investment costs are much lower compared to, for example, optical measurement systems. The wide range of different sieving methods is fully covered by RETSCH's offering which provides a suitable sieve shaker for practically any bulk material. RETSCH sieve shakers guarantee exact and reproducible results in a very short time and fulfill all requirements for test agent monitoring according to DIN EN ISO 9000 ff.

Expert Guide Sieving

Would you like to learn more about sieve analysis?
Then download the 50-page expert guide from our website!

www.retsch.com/downloads/brochures

